

INSTRUMENTAL ANALYSIS II



Unit 4

UV-Vis Spectroscopy


Course Instructor: Ermias Haile (Ass.Prof.)

INTRODUCTION

- UV-VIS spectroscopy is the measurement of the wavelength and intensity of absorption of near-ultraviolet and visible light by a sample.
- Ultraviolet and visible light are energetic enough to promote outer electrons to higher energy levels.
- UV-VIS spectroscopy is usually applied to molecules and inorganic ions or complexes in solution.
- UV-VIS spectra have broad features that are of limited use for sample identification but are very useful for quantitative measurements. The concentration of an analyte in solution can be determined by measuring the absorbance at some wavelength and applying the Beer-Lambert Law

UV radiation and Electronic Excitations


- ✓ The **difference in energy** between **molecular bonding**, **non-bonding** and **anti-bonding** orbitals ranges from **125-650 kJ/mole**
- ✓ This energy corresponds to EM radiation in the **ultraviolet (UV) region, 100-350 nm**, and **visible (VIS) regions 350-750 nm** of the spectrum
- ✓ For comparison, recall the EM spectrum:



- ✓ Using **IR** will observe **vibrational** transitions with energies of **8-40 kJ/mol** at wavelengths of **2500-15,000 nm**
- ✓ For purposes of our discussion, we will refer to UV and VIS spectroscopy as UV

The Spectroscopic Process

- In UV spectroscopy, the sample is irradiated with the broad spectrum of the UV radiation
- If a particular electronic transition matches the energy of a certain band of UV, it will be absorbed
- The remaining UV light passes through the sample and is observed
- From this residual radiation a spectrum is obtained with “gaps” at these discrete energies-this is called an absorption spectrum




Observed electronic transitions

- ❖ The lowest energy transition most often observed by UV is typically that of an electron in the Highest Occupied Molecular Orbital (HOMO) to the Lowest Unoccupied Molecular Orbital (LUMO)
- For any bond (pair of electrons) in a molecule, the molecular orbitals are a mixture of the two contributing atomic orbitals; for every bonding orbital “created” from this mixing (σ , π), there is a corresponding anti-bonding orbital of symmetrically higher energy (σ^* , π^*)
- ❖ The lowest energy occupied orbitals are typically the σ ; likewise, the corresponding anti-bonding σ^* orbital is of the highest energy
- ❖ π -orbitals are of somewhat higher energy, and their complementary anti-bonding orbital somewhat lower in energy than σ^* .
- ❖ Unshared pairs lie at the energy of the original atomic orbital, most often this energy is higher than π or σ .

Cont...

□ Here is a graphical representation



□ From the molecular orbital diagram, there are several possible electronic transitions that can occur, each of a different relative energy:

$\sigma \rightarrow \sigma^*$ alkanes


\rightarrow

$\pi \rightarrow \pi^*$ unsaturated cpds.

$n \rightarrow \sigma^*$ O, N, S, halogens

$n \rightarrow \pi^*$ carbonyls


An Electronic Spectrum



By Ermias H.


Absorption Spectrum.

An **absorption spectrum** is the plot of the **absorbance** versus the **wavelength** of the incident light.



Transmission Spectrum.

A transmittance spectrum is the plot of the transmittance or %T versus the wavelength of the incident light.



Types of Electronic Transitions

1. Transitions involving π , σ , and n electrons
2. Transitions involving charge-transfer electrons
3. Transitions involving d and f electrons (not covered in this Unit)

❖ **Four** types of transitions are possible:

$\sigma \rightarrow \sigma^*$, $n \rightarrow \sigma^*$, $n \rightarrow \pi^*$, and $\pi \rightarrow \pi^*$.




Figure : Electronic molecular energy levels

Absorbing species containing π , σ , and n electrons

- Absorption of ultraviolet and visible radiation in organic molecules is restricted to certain functional groups (*chromophores*) that contain valence electrons of low excitation energy
- The spectrum of a molecule containing these chromophores is complex.
- This is because the superposition of rotational and vibrational transitions on the electronic transitions gives a combination of overlapping lines.
- This appears as a continuous absorption band.

$\sigma - \sigma^*$ Transitions

- An electron in a bonding σ orbital is excited to the corresponding antibonding orbital.
- The energy required is large.
- For example, methane (which has only C-H bonds, and can only undergo $\sigma - \sigma^*$ transitions) shows an absorbance maximum at 125 nm.
- Absorption maxima due to $\sigma - \sigma^*$ transitions are not seen in typical UV-Vis. spectra (200 - 700 nm)

$n - \sigma^*$ Transitions


- Saturated compounds containing atoms with lone pairs (non-bonding electrons) are capable of $n - \sigma^*$ transitions.
- These transitions usually need less energy than $\sigma - \sigma^*$ transitions. They can be initiated by light whose wavelength is in the range 150 - 250 nm.
- The number of organic functional groups with $n - \sigma^*$ peaks in the UV region is small.

$n - \pi^*$ and $\pi - \pi^*$ Transitions

- Most absorption spectroscopy of organic compounds is based on transitions of n or π electrons to the π^* excited state.
- This is because the absorption peaks for these transitions fall in an experimentally convenient region of the spectrum (200 - 700 nm).
- These transitions need an unsaturated group in the molecule to provide the π electrons.

Vacuum UV or Far UV ($\lambda < 190 \text{ nm}$)


UV/VIS



Example


❖ **Ethanal** can therefore absorb light of two different wavelengths:

- the π bonding to π^* anti-bonding absorption peaks at 180 nm;
- the n non-bonding to π^* anti-bonding absorption peaks at 290 nm.



The importance of conjugation and delocalisation in what wavelength is absorbed

❖ Consider these **three** molecules:



ethene

buta-1,3-diene


hexa-1,3,5-triene

- Ethene contains a simple isolated carbon-carbon double bond, but the other two have conjugated double bonds.
- In these cases, there is delocalisation of the pi bonding orbitals over the whole molecule.
- The maximum absorption is moving to **longer wavelengths as the amount of delocalisation increases.**


molecule	wavelength of maximum absorption (nm)
ethene	171
buta-1,3-diene	217
hexa-1,3,5-triene	258

❖ Compare ethene with buta-1,3-diene.

- In **ethene**, there is one pi bonding orbital and one pi anti-bonding orbital.
- In **buta-1,3-diene**, there are two pi bonding orbitals and two pi anti-bonding orbitals.
- The highest occupied molecular orbital is often referred to as the **HOMO** – in these cases, it is a pi bonding orbital.
- The lowest unoccupied molecular orbital (the **LUMO**) is a pi anti-bonding orbital.



❖ In the hexa-1,3,5-triene case, it is **less** still.



➤ If you extend this to compounds with really massive delocalization, the wavelength absorbed will eventually be high enough to be in the **visible region** of the spectrum, and the compound will then be seen as coloured.

The UV/visible spectrum

- Figure shows a simple UV-visible spectrum.
- Absorbance (on the vertical axis) is a measure of the amount of light absorbed.
- The higher the value, the more of a particular wavelength is being absorbed.





Figure : Typical absorption spectra for buta-1,3-diene ($\text{CH}_2=\text{CH}-\text{CH}=\text{CH}_2$)

Typical Absorptions of Simple Isolated Chromophores

Class	Transition	λ_{\max} (nm)	$\log \epsilon$	Class	Transition	λ_{\max} (nm)	$\log \epsilon$
R—OH	$n \rightarrow \sigma^*$	180	2.5	R—NO ₂	$n \rightarrow \pi^*$	271	<1.0
R—O—R	$n \rightarrow \sigma^*$	180	3.5	R—CHO	$\pi \rightarrow \pi^*$	190	2.0
R—NH ₂	$n \rightarrow \sigma^*$	190	3.5		$n \rightarrow \pi^*$	290	1.0
R—SH	$n \rightarrow \sigma^*$	210	3.0	R ₂ CO	$\pi \rightarrow \pi^*$	180	3.0
R ₂ C=CR ₂	$\pi \rightarrow \pi^*$	175	3.0		$n \rightarrow \pi^*$	280	1.5
R—C≡C—R	$\pi \rightarrow \pi^*$	170	3.0	RCOOH	$n \rightarrow \pi^*$	205	1.5
R—C≡N	$n \rightarrow \pi^*$	160	<1.0	RCOOR'	$n \rightarrow \pi^*$	205	1.5
R—N=N—R	$n \rightarrow \pi^*$	340	<1.0	RCONH ₂	$n \rightarrow \pi^*$	210	1.5

Instrumentation and Spectra

- ❖ The construction of a traditional **UV-VIS spectrometer** consists — sample handling, irradiation, detection and output
- Here is a simple schematic that covers **most modern UV spectrometers**



- ❖ **Two sources** are required to scan the entire UV-VIS band:
 - Deuterium** lamp — covers the UV — 200-330
 - Tungsten** lamp — covers 330-700



Sample Handling

- almost all UV spectra are recorded **solution-phase**
- Cells can be **made** of plastic, **glass** or quartz
- Only **quartz** is transparent in the full **200-700 nm range**;
- **plastic** and **glass** are only suitable for **visible spectra**
- Concentration (we will cover shortly) is empirically determined
- A typical sample cell (commonly called a **cuvet**):




Solvents

- ❖ The choice of the solvent to be used in UV spec. is quite important.
- The first criterion for a good solvent is that it **should not be absorb UV radiation**(must be transparent) in the same region as the substance whose spectrum is being determined.
- Usually solvents that **do not contain conjugated systems** are most suitable for this purpose.
- In the **UV**, any solvent that does not have any π bonding can be used including water.
- In the **visible** region any solvent that is colorless can be used which also includes water.
- Each is transparent in the regions of UV spectrum where interesting absorption peaks from sample molecules are likely to occur.

The Spectrum

- The **x-axis** of the spectrum is in **wavelength**; 200-350 nm for UV, **200-700** for UV-VIS determinations, The **y-axis** of the spectrum is in **absorbance, A**
- Due to the lack of any fine structure, spectra are rarely shown in their raw form, rather, the peak maxima are simply reported as a numerical list of “lambda max” values or λ_{max}



$\lambda_{\text{max}} =$ 206 nm
252
317
376

Cont...

- ❖ From the spectrometers point of view, **absorbance is the inverse of transmittance**: $A = \log_{10} (I_0/I)$
- ❖ From an **experimental point of view**, **three** other **considerations** must be made:
 - i. A **longer path length**, l through the sample will **cause more UV light to be absorbed** – linear effect
 - ii. The **greater** the **concentration**, c of the sample, the **more UV light will be absorbed** – linear effect
 - iii. Some **electronic transitions are more effective** at the absorption of photon than others – **molar absorptivity**, ϵ this may vary by orders of magnitude.
- ❖ These effects are combined into the Beer-Lambert Law: $A = \epsilon c l$
- ❖ For **most UV** spectrometers, l would remain **constant** (standard cells are typically **1 cm in path length**)

Practical application of UV spectroscopy

- UV spectra and visible spectra can be used to identify an **unknown compound** by a comparative analysis
 - UV was the first organic spectral method, however, it is rarely used as a primary method **for structure determination**
 - It is most **useful in combination with NMR and IR data** to elucidate unique electronic features that may be ambiguous in those methods
 - It can be used to assay (via λ_{max} and molar absorptivity) the proper irradiation λ s for photochemical experiments, or the design of UV resistant paints and coatings
 - The most everywhere use of **UV is as a detection device for HPLC**; since UV is utilized for solution phase samples vs. a reference solvent this is easily incorporated into LC design
- * UV is detector to HPLC as mass spectrometry (MS) to GC

Chromophores


A. Definition

- Remember the electrons present in organic molecules are involved in covalent bonds or lone pairs of electrons on atoms such as O/N
- Since similar functional groups will have es capable of discrete classes of transitions, the characteristic energy of these energies is more representative of the functional group than the es themselves
- A functional group capable of having characteristic electronic transitions is called a chromophore (color loving)
- Structural or electronic changes in the chromophore can be quantified and used to predict shifts in the observed electronic transitions

B. Organic Chromophores


1. **Alkanes** – only possess σ -bonds and no lone pairs of electrons, so only the high energy $\sigma \rightarrow \sigma^*$ transition is observed in the far UV

- This transition is destructive to the molecule, causing cleavage of the σ -bond



2. Alcohols, ethers, amines & sulfur cpds


- In the cases of simple, aliphatic examples of these cpds the $n \rightarrow \sigma^*$ is the most often observed transition; like the alkane $\sigma \rightarrow \sigma^*$ it is most often at shorter λ than 200 nm. How this transition occurs from the HOMO to the LUMO



27

3. Alkenes & Alkynes

- In the case of isolated examples of these cpds the $\pi \rightarrow \pi^*$ is observed at 175 and 170 nm, respectively
- Even though this transition is of lower energy than $\sigma \rightarrow \sigma^*$ it is still in the far UV- however, the transition energy is sensitive to substitution




4. Carbonyls

- Unsaturated systems incorporating N or O can undergo $n \rightarrow \pi^*$ transitions (~ 285 nm) in addition to $\pi \rightarrow \pi^*$
- Despite the fact this transition is forbidden by the selection rules ($\epsilon = 15$), it is the most often observed & studied transition for carbonyls
- This transition is also sensitive to substituents on the carbonyl
- Similar to alkenes and alkynes, non-substituted carbonyls undergo the $\pi \rightarrow \pi^*$ transition in the vacuum UV (188 nm, $\epsilon = 900$); sensitive to substitution effects

Cont...

Carbonyls – $n \rightarrow \pi^*$ transitions (~ 285 nm); $\pi \rightarrow \pi^*$ (188 nm)




- It has been determined from spectral studies, that carbonyl oxygen more approximates sp rather than sp^2 !
- σ_{CO} transitions omitted for clarity

C. Substituent Effects







- ❖ **General** – from our brief study of these general chromophores, only the weak $n \rightarrow \pi^*$ transition occurs in the routinely observed UV
- The attachment of **substituent** groups (other than H) **can shift the energy of the transition**
- **Substituents** that **increase** the **intensity** and often **wavelength** of an absorption are called **auxochromes**
- **Common auxochromes** include **alkyl**, **hydroxyl**, **alkoxy** and **amino** groups and the **halogens**

- ❖ Generally- **Substituents** may have any of **four** effects on a chromophore
- Bathochromic shift** (red shift) – a **shift to longer λ** ; lower energy
 - Hypsochromic shift** (blue shift) – **shift to shorter λ** ; higher energy
 - Hyperchromic effect** – an **increase in intensity**
 - Hypochromic effect** – a **decrease in intensity**




Cont...

1. Conjugation- most efficient means of bringing about a bathochromic (shift to longer λ) and hyperchromic shift of an unsaturated chromophore:

	λ_{max} nm	ϵ
	175	15,000
	217	21,000
	258	35,000
 β -carotene	465	125,000
	$n \rightarrow \pi^*$ 280	12
	$\pi \rightarrow \pi^*$ 189	900
	$n \rightarrow \pi^*$ 280	27
	$\pi \rightarrow \pi^*$ 213	7,100


Cont...

- ❖ **Conjugation – Alkenes** : The observed shifts from conjugation imply that an increase in conjugation decreases the energy required for electronic excitation
- ❖ Extending this effect out to longer conjugated systems the energy gap becomes progressively smaller:




Colour and light absorption- The Chromophore concept

- ❖ Compounds that absorb light of wavelength between 400 and 800 nm (visible light) appear colored to the human eye.
- ❖ When white light passes through a solution and come into view as red light, we say that the solution is red.
- ❖ What has actually happened is that the solution has allowed the red component of white light to pass through, whereas it has absorbed the complementary colors, yellow and blue.
- ❖ The more concentrated the sample solution, the more yellow and blue light is absorbed and the more intensely red the solution appears to the eye.



The human eye sees the complementary color to that which is absorbed

Absorbance and Complementary Colors



Visible Spectroscopy

Color

General:- The portion of the EMS from 400-800nm is observable to humans. We & some other mammals have the adaptation of seeing color at the expense of greater detail




Violet	400-420
Indigo	420-440
Blue	440-490
Green	490-570
Yellow	570-585
Orange	585-620
Red	620-780

V. Visible Spectroscopy

A. Color

1. General

- When **white** (continuum of λ) light passes through, or is reflected by a surface, those that are absorbed are removed from the transmitted or reflected light respectively
- What is “**seen**” is the **complimentary colors** (those that are not absorbed)
- This is the origin of the “color wheel”





V. Visible Spectroscopy

A. Color

1. General

- Organic compounds that are “colored” are typically those with extensively **conjugated** systems (typically **more** than **five**)

- Consider β -carotene








λ_{max} is at 455 – in the far **blue** region of the spectrum – this is **absorbed**

The remaining light has the **complementary** color of **orange**

Substituent Effects

1. **Conjugation** – most efficient means of bringing about a bathochromic and hyperchromic shift of an unsaturated chromophore:

$\text{H}_2\text{C}=\text{CH}_2$		λ_{max} nm	ϵ
		175	15,000
		217	21,000
		258	35,000
 β -carotene		465	125,000
	$n \rightarrow \pi^*$	280	12
	$\pi \rightarrow \pi^*$	189	900
	$n \rightarrow \pi^*$	280	27
	$\pi \rightarrow \pi^*$	213	7,100

IV. Structure Determination

A. Dienes

2. Woodward-Fieser Rules

Woodward and the Fiesers performed extensive studies of terpene and steroidal alkenes and noted similar substituents and structural features would predictably lead to an empirical prediction of the wavelength for the lowest energy $\pi \rightarrow \pi^*$ electronic transition

This work was distilled by Scott in 1964 into an extensive treatise on the Woodward-Fieser rules in combination with comprehensive tables and examples – (A.I. Scott, Interpretation of the Ultraviolet Spectra of Natural Products, Pergamon, NY, 1964)

A more modern interpretation was compiled by Rao in 1975 – (C.N.R. Rao, Ultraviolet and Visible Spectroscopy, 3rd Ed., Butterworths, London, 1975)

❖ This quantification is referred to as the **Woodward-Fieser** Rules which we will apply to **three specific chromophores**:

1. **Conjugated dienes**
2. **Conjugated dienones**
3. **Aromatic systems**


A. Dienes

1. General Features

For acyclic butadiene, two conformers are possible – s-cis and s-trans




s-trans



s-cis


- The s-cis conformer is at an overall higher potential energy than the s-trans; therefore the HOMO electrons of the conjugated system have less of a jump to the LUMO – lower energy, longer wavelength




$\lambda_{\max} = 217$




253




220




227



227



256



263 nm

A. Dienes

2. Woodward-Fieser Rules - Dienes

The rules begin with a base value for λ_{max} of the chromophore being observed:



acyclic butadiene = 217 nm


The incremental contribution of substituents is added to this base value from the group tables:

Group	Increment
Extended conjugation	+30
Carbo cyclic ring	+5
Alkyl	+5
-OCOCH ₃	+0
-OR	+6
-SR	+30
-Cl, -Br	+5
-NR ₂	+60

A. Dienes

2. Woodward-Fieser Rules - Dienes


For example:



Isoprene - acyclic butadiene =
one alkyl subs.

Experimental value

217 nm
+ 5 nm
222 nm
220 nm



Allylidencyclohexane
- acyclic butadiene =
one exocyclic C=C
2 alkyl subs.

Experimental value


217 nm
+ 5 nm
+10 nm
232 nm
237 nm

A. Dienes

3. Woodward-Fieser Rules – Cyclic Dienes


There are two major types of cyclic dienes, with two different base values

Heteroannular (transoid):



$$\varepsilon = 5,000 - 15,000$$
$$\text{base } \lambda_{\text{max}} = 214$$

Homoannular (cisoid):



$$\varepsilon = 12,000 - 28,000$$
$$\text{base } \lambda_{\text{max}} = 253$$


The increment table is the same as for acyclic butadienes with a couple additions:

Group	Increment
Additional homoannular	+39
Where both types of diene are present, the one with the longer λ becomes the base	


A. Dienes

3. Woodward-Fieser Rules – Cyclic Dienes


For example:




1,2,3,7,8,8a-hexahydro-8a-methylnaphthalene



heteroannular diene = 214 nm



3 alkyl subs. (3 x 5) +15 nm




1 exo C=C

+ 5 nm
234 nm

Experimental value 235 nm

A. Dienes


3. Woodward-Fieser Rules – Cyclic Dienes



heteroannular diene = 214 nm

4 alkyl subs. (4 x 5) +20 nm
1 exo C=C + 5 nm

239 nm



homoannular diene = 253 nm


4 alkyl subs. (4 x 5) +20 nm
1 exo C=C + 5 nm

278 nm


A. Dienes

3. Woodward-Fieser Rules – Cyclic Dienes


Be careful with your assignments – three common **errors**:



This compound has three exocyclic double bonds; the indicated bond is exocyclic to two rings



This is **not a heteroannular** diene; you would use the base value for an acyclic diene



Likewise, this is **not a homoannular** diene; you would use the base value for an acyclic diene